

Simple Test Process

1 Scan or enter operator and patient ID



2 Insert strip into meter



3 Prick finger to obtain blood



4 Apply blood drop



5 Read results in 30 seconds



Coming Soon!

Connectivity to hospital network or **IRISCT™**

ORDERING INFORMATION

Description	Cat. No.
EZ Chem Meter Kit	9600K
EZ Chem Meter	9600
Docking Station with AC Adaptor	9602
EZ Chem Creatinine Test Strips, 25 per vial, 2 vials per package	9604
Battery for EZ Chem Meter, 4-Pack	9606
Carrying Case for Meter & Supplies	9608
EZ Chem Creatinine Control Solution, 3 Level Pack	9612
EZ Chem Creatinine Linearity Kit, 5 levels, 1 vial of each level	9618
Instructions for use manual, printed	9620
Instructions for use manual, CD	9622
Quick Reference Guide	9624

To order additional supplies, contact Customer Solutions at 800-544-4624.



For information regarding the acquisition of E-Z-EM by Bracco Diagnostics Inc., visit our homepage at www.ezem.com.